



## Patient Simulator

*Psi 300*



## Instruction Manual



© 04/2007, Peter Kron GmbH  
Subject to Technical Changes

**mtk Peter Kron GmbH Mess-Systeme für Medizintechnik Hard & Software**  
Zossener Straße 41, D-10961 Berlin, Germany  
Tel.: +49 (0)30 / 69 81 88 - 40 Fax: +49 (0)30 69 81 88 - 49

## Contents

|   |    |
|---|----|
| 1. Purpose and Application .....                        | 4  |
| 2. Models.....  | 4  |
| 3. Delivery .....                                       | 5  |
| 4. Optional Accessories .....                           | 5  |
| 5. Important Security Notice.....                       | 6  |
| 6. Power Supply .....                                   | 6  |
| 7. Operating Keys .....                                 | 8  |
| 8. Display Elements .....                               | 9  |
| 9. Connections .....                                    | 10 |
| 10. Getting Started .....                               | 11 |
| 11. Operating Instructions .....                        | 12 |
| 12. Operational Sequences and Definitions.....          | 16 |
| 12.1 Main Menu .....                                    | 16 |
| 12.2 Main Function "ECG" .....                          | 17 |
| 12.3 Main Function "RESP" - Respiration.....            | 21 |
| 12.4 Main Function "BP" – Invasive Blood Pressure ..... | 22 |
| 12.5 Main Function "TEMP" - Temperature.....            | 24 |
| 12.6 Main Function "C.O." – Cardiac Output .....        | 26 |
| 12.7 Main Function "NS" - Noise.....                    | 27 |
| 13. Environmental Factors .....                         | 28 |
| 14. Cleaning and Disinfection .....                     | 28 |
| 15. Errors, Defects and Customer Service .....          | 29 |
| 15.1 Error Codes.....                                   | 29 |
| 15.2 Trouble-Shooting .....                             | 29 |

|   |    |
|---|----|
| 16. Maintenance.....                                      | 29 |
| 16.1 Updating the Operating Software (Firmware Update)... | 30 |
| 16.2 Determining the Current Software .....               | 30 |
| 17. Intervals for Calibration .....                       | 30 |
| 18. Packaging Materials and Disposal.....                 | 30 |
| 19. Pictographs .....                                     | 31 |
| 20. Technical Specifications.....                         | 32 |

## 1. Purpose and Application

The PSI 300 is a patient simulator. It generates signals that simulate predetermined sequences in a corresponding recording apparatus, such as a patient monitor, and thereby to interpret and display the patient's vital parameter.

The PSI 300 simulates the following vital parameter:

- Electrocardiogram           **[ ECG ]**
- Respiration                   **[ RESP ]**
- Invasive blood pressure   **[ BP ]**
- Temperature                 **[ TEMP ]**
- Cardiac output               **[ C.O. ]** *with PSI 300 only – Extended version*

## 2. Models

The PSI 300 is available in 4 models that differ in function according to the respective software (firmware).

### • PSI 300 – Basic Model

The basic model of the PSI 300 with the firmware "V 1.xx" comprises the simulation of 4 vital parameter:

- *Electrocardiogram, 12 channels*
- *Respiration*
- *Invasive blood pressure*
- *Temperature*

### • PSI 300 – Basic Model With RS-232 Interface

The basic model of the PSI 300 with the firmware version "RS 1.xx" offers the above-mentioned 4 basic vital parameter simulation, but also includes the RS-232 interface for connection to a PC.

### • PSI 300 – Extended Model "C.O."

The extended version of the PSI 300 with firmware "V 1.xx-CO" allows the simulation of 5 vital parameter:

- *Electrocardiogram, 12 channels*
- *Respiration*

- *Invasive blood pressure*
- *Temperature*
- *Cardiac output*

### • **PSI 300 – Extended Model "C.O." With RS-232 Interface**

The extended model of the PSI 300 with the firmware version "RS 1.xx-CO" offers not only the above-mentioned 5 basic vital parameter simulation, but also the RS-232 interface for connection to a PC as well.

#### **Note:**

After turning on the PSI 300, during the initialization, the respective version and firmware revision (xx) numbers are shown in the display for approximately 2 seconds.

### **3. Delivery**

The shipment of the PSI 300 includes the following:

- 1 Patient simulator PSI 300, Item No. 8210 0014
- 4 Batteries 1.5V (AA)
- 1 Nullmodem-cable RS 232  
(only with firmware versions "RS 1.xx" / "RS 1.xx-CO")
- 1 Instruction manual
- 1 Manufacturer's Factory Calibration Certificate  
*Please deposit this certificate in a safe place, together with your future calibration protocols.*

### **4. Optional Accessories**

The following optional accessories are available for the PSI 300

| Item            | Description  |
|-----------------|--|
| Battery charger | Charger for NiCd or NiMH rechargeable batteries (15 minutes) |

## 5. Important Security Notice

Please read the instructions in this manual carefully before using the Patient Simulator PSI 300 the first time. The manufacturer does not assume liability or a warranty for mistakes or damages resulting from incorrect operation or use.

- Never use the PSI 300 directly on a patient or on an instrument that is currently being used with a patient.
- Place the PSI 300 so that the LCD display is clearly visible. Avoid direct sunlight.
- Do not press on the LCD display.
- Avoid short circuiting the signals on the output jacks. Never use damaged cables.
- Make sure no liquids penetrate the inside of the PSI 300.
- In case of a defect, remove the batteries and contact your Customer Service. Only trained and qualified personnel are authorized to repair the PSI 300. If an unauthorized person opens the apparatus, the warranty and any liability claims automatically expire.
- Keep this Instruction Manual handy for future reference.

## 6. Power Supply

The PSI 300 can be operated with commercial alkaline batteries or rechargeable storage batteries.

### • Batteries

4 Batteries á 1.5V (AA or LR6), Alkaline

The average operating time with one set of batteries is ca. 70 hours.

### • Rechargeable Storage Batteries

4 Storage Batteries á 1.2V (AA or LR6), NiCd or NiMH

The average operating time with one set of NiMH, 1800mAh storage batteries is ca. 50 hours.

The storage batteries must be recharged in an appropriate external charger. A charger with a loading time of only 15 minutes is available as an optional accessory.

## Caution!

Change the batteries as described in Section 10, if the any of the following occurs:

the LCD display has become weak and doesn't have a clear contrast,  
the PSI 300 shuts off automatically,  
the PSI 300 cannot be turned on.

*In these cases, immediately remove the weak batteries!*

**Note:**



Never dispose of used batteries and storage batteries in ordinary household trash. Adhere to your local regulations and laws for disposal. Always recycle used batteries.

## 7. Operating Keys

All operating elements of the PSI 300 are located on the keyboard:



"ON"



"OFF"



Cursor key "UP"



Cursor key "DOWN"



Cursor key "LEFT"



Cursor key "RIGHT"



"EXECUTE" the setting using the cursor keys



Function key for choosing main functions:






## 8. Display Elements


The current operating conditions and all settings are shown in the four-line **LCD-Monitor** à 40 characters.

Using the cursor keys, the contents of the window can be shifted up or down as well as to the left or right, according to the main function.



## 9. Connections

The DIN jacks for connecting the cables for:  
invasive blood pressure,  
temperature, and  
cardiac output  
are located on the back of the PSI 300.



The snap fasteners for connecting the  
ECG signal and  
respiration  
are located on the upper part of the front of the PSI 300.

### \*Note:

Cardiac output can only be simulated with the:  
*PSI 300 – Extended version "C.O."* or  
*PSI 300 – Extended version "C.O." with RS-232 interface*  
together with the firmware "V 1.xx-CO" or "RS 1.xx-CO," respectively.


### \*\*Note:

The Interface RS-232 for connecting the PSI 300 to a PC is supported only  
with the firmware "RS 1.xx" or "RS 1.xx-CO," respectively.

## 10. Getting Started

Before you start, insert new batteries or fully charged storage batteries. See the instructions in Section 6.

- In order to insert the batteries, open the cover of the battery compartment on the bottom of the PSI 300.




Put the tip of a small screwdriver into the smaller slot of the cover and gently push the cover up. The cover will open easily.

### Caution!

Pay strict attention to the right position of the batteries (poles)!

- **Close the cover of the battery compartment after inserting the batteries.**



Insert the two pegs of the cover into the holes of the battery compartment and gently push the cover down until it clicks into place.

### Caution!

The cover can only be inserted in one direction.

- **Turn on the PSI 300 as described in Section 11 below.**

## 11. Operating Instructions

All operations of the PSI 300 are managed by menus and are executed using 13 menu keys on the keyboard.

Operating the PSI 300 via the interface is only possible with those models that run under the firmware "RS 1.xx" or "RS 1.xx-CO."

- **Starting the PSI 300**

Switch on the PSI 300 by pressing the **"On" button**.

After starting, during the initialization process, the model and the firmware version of the PSI 300 are shown in the display for about 2seconds.

|                                   |
|-----------------------------------|
| PSI 300 Patient Simulator V1.0-CO |
| power-on selftest running .....   |

When the Start menu appears in the display, the PSI 300 is ready for use.

|                               |
|-------------------------------|
| [ECG] RESP BP TEMP C.O. NS    |
| ECG form [IEC 601] others     |
| type [calibration] analytical |
| heart rate [60] 120           |

**Note:**

Every time the PSI 300 is switched on, it will start with the main functions default. The next menu to appear is the Start menu for the main function "ECG."

- **Turning Off the PSI 300**

Turn off the PSI 300 by pressing the **"Off" key**. Switch i off when you are not going to use it any longer period of time.

**Note:**

Any changes in the main functions are not stored and need to be set up again each time you switch on the PSI 300.

## • Choosing the Main Function

In order to change to another function shown in the Main Menu (display Line 1), press the **Function Key** that is located above the respective function.



|            |               |            |      |      |    |
|------------|---------------|------------|------|------|----|
| [ECG]      | RESP          | BP         | TEMP | C.O. | NS |
| ECG form   | [IEC 601]     | others     |      |      |    |
| type       | [calibration] | analytical |      |      |    |
| heart rate | [60]          | 120        |      |      |    |

Line 1

Main function chosen in this example: **[ECG]**

### Note:

The main menu in Line 1 always appears in the monitor.

## Subfunctions and Parameter

The subfunctions and parameter of the selected main function are shown in Lines 2 to 4.

|            |               |            |      |      |    |
|------------|---------------|------------|------|------|----|
| [ECG]      | RESP          | BP         | TEMP | C.O. | NS |
| ECG form   | [IEC 601]     | others     |      |      |    |
| type       | [calibration] | analytical |      |      |    |
| heart rate | [60]          | 120        |      |      |    |

Line 1

Line 2

Line 3

Line 4

Active subfunction in this example: **[artificial], [calibration]**

Active parameter: **[60]**

## Cursor Key s"Down" and "Up"



For the main functions "ECG" and "RESP," further lines can be shown by pressing the cursor key "Down." The display shifts down one line each time the key is pressed:

|            |               |            |      |         |    |
|------------|---------------|------------|------|---------|----|
| [ECG]      | RESP          | BP         | TEMP | C.O.    | NS |
| type       | [calibration] | analytical |      |         |    |
| heart rate | [60]          | 120        |      |         |    |
| QRS form   | [RS]          | R          | QS   | smallRS |    |

Line 1

Line 3

Line 4

Line 5

Press the cursor key "Up" in order to return to a line further up the list. The display will move up one line each time the key is pressed, until Lines 2 to 4 reappear under Line 1.

**Note: The Main Menu in Line 1 is always shown in the display.**

## Cursor Keys "Right" and "Left"



For the main functions "TEMP" and "C.O.," further parameter can be selected by pressing the cursor key "Right" if the symbol ">" appears at the end of the line. The display shifts one field further to the right each time the cursor key "Right" is pressed.

|            |         |       |       |        |         |
|------------|---------|-------|-------|--------|---------|
| ECG        | RESP    | BP    | TEMP  | [C.O.] | NS      |
| Blood temp | (°C)    | 36    | [37]  | 38     |         |
| min vol    | (l/min) | < 4.5 | [5.0] | 5.5    | 6.0 7.> |

Field 41...

Press the cursor key "Left" in order to return to the previous parameter fields that are no longer visible. The display will shift back one field further to the left each time the cursor key "Left" is pressed, until the first parameter field is back in the first position.

### Selecting the Subfunctions and Changing the Parameter

The active subfunctions and parameter are shown in brackets.


|            |               |            |      |      |    |
|------------|---------------|------------|------|------|----|
| [ECG]      | RESP          | BP         | TEMP | C.O. | NS |
| ECG form   | [IEC 601]     | others     |      |      |    |
| type       | [calibration] | analytical |      |      |    |
| heart rate | [60]          | 120        |      |      |    |

Line 1  
 Line 2  
 Line 3  
 Line 4



Active subfunctions in this example: **[artificial], [calibration]**

Active parameter: **[60]**



Using the yellow **direction keys** (up, down, left, right), move the cursor in the display field over the first letter of the desired subfunction or the parameter.

Complete your choice by pressing the green **"Execute"** key.

### CAUTION!



In the main functions "ECG" and "RESP," more lines are available in order to change the functions or the parameter by pressing the "Down" cursor key.

The display moves one line lower each time the "Down" cursor key is pressed. The main menu in Line 1 always remains in view.

Press the "Up" cursor key in order to return to the upper lines. The display moves one line higher each time the "Up" cursor key is pressed, until Lines 1, 2, 3, and 4 are shown in the display.

**Note:**

After turning on the PSI 300, or after making changes, the functions and parameter are still active when a different menu appears on the display after pressing a function selection key. They remain active until a new function is implemented by pressing the “Execute” key.

## 12. Operational Sequences and Definitions

### Main Menu

The main menu is always shown in the **first (top) line** of the display when the PSI 300 is in operation.

|              |               |            |             |             |           |
|--------------|---------------|------------|-------------|-------------|-----------|
| <b>[ECG]</b> | <b>RESP</b>   | <b>BP</b>  | <b>TEMP</b> | <b>C.O.</b> | <b>NS</b> |
| ECG form     | [IEC 601]     | others     |             |             |           |
| type         | [calibration] | analytical |             |             |           |
| heart rate   | [60]          | 120        |             |             |           |

The meanings and explanations of the abbreviations shown in the main menu are the following:

| Display      | Meaning           | Selection (Function Key) |
|--------------|-------------------|--------------------------|
| <b>[ECG]</b> | Electrocardiogram | Main function "ECG"      |
| <b>RESP</b>  | Respiration       | Main function "RESP"     |
| <b>BP</b>    | Blood Pressure    | Main function "BP"       |
| <b>TEMP</b>  | Temperature       | Main function "TEMP"     |
| <b>C.O.</b>  | Cardiac Output    | Main function "C.O."     |
| <b>NS</b>    | Noise             | Main function "NS"       |



## Main Function "ECG"

In the default setting after turning on the PSI 300, the main function "ECG" appears in the top line of the display, together with the following settings in Lines 2, 3 and 4:

|            |               |            |      |      |    |
|------------|---------------|------------|------|------|----|
| [ECG]      | RESP          | BP         | TEMP | C.O. | NS |
| ECG form   | [IEC 601]     | others     |      |      |    |
| type       | [calibration] | analytical |      |      |    |
| heart rate | [60]          | 120        |      |      |    |

Line 1 (Main Menu)

Line 2

Line 3

Line 4



### Note:

The settings and parameter in Line 5 to Line 8 can be shown and changed by pressing the "down" cursor key.

The basic settings of the main function "ECG" can be changed according to the following table:

| Line 2                    | Line 3              | Line 4  | Line 5  | Line 6   | Line 7  | Line 8              |
|---------------------------|---------------------|---|---|--|---|---------------------|
| ECG form:<br>[IEC 601]    | type:<br>[calibrat] | heart rate<br>[60]                                | QRS form:<br>[RS]   | QRS(mV):<br>±0.5 ±1.0 ±1.5<br>[±2.0] ±3.0 ±5.0                       | STdev(µV):<br>[0]   | -                   |
|                           |                     |   | QRS form:<br>[R] QS   | QRS(mV):<br>[±2.0]   | STdev(µV):<br>-200 [0] +200   | -                   |
|                           |                     |   | QRS form:<br>[smallRS]  | QRS(mV):<br>[±2.0]   | STdev(µV):<br>[0]   | -                   |
|                           | type:<br>analytical | heart rate<br>[120]                               | QRS form:<br>[RS]   | QRS(mV):<br>[±2.0]   | STdev(µV):<br>[0]   | -                   |
| heart rate<br>40 [60] 120 |                     | QRS form:<br>[RS]                                 | QRS(mV):<br>±2.0  | STdev(µV):<br>[0]  | -   |                     |
| ECG form:<br>others       | [sinusrh]           | heart rate<br>30 [60]<br>80 120<br>160 200<br>240 | Axis:<br>[INT] HOR VERT   | QRS(mV):<br>0.05 0.1 0.2 0.5<br>0.7 [1.0] 1.5 2.0<br>2.5 3.0 4.0 5.0 | STdev(µV):<br>-800 -600 -400 -<br>200 -100 -50 [0]<br>+50 +100 +200<br>+400 +600 +800 | NEO:<br>[OFF]<br>ON |
|                           |                     | arrhyth   | category<br>[suprve]  | type<br>[af1] af2 af1t sa<br>mb80 mb120<br>pat nodal svt             | -   | -                   |
|                           | category<br>[prem]  |   | type<br>[pac] pnc pvc<br>mf1 mf2 mf3                                | -  | -   | -                   |
|                           | category<br>[vntar] |   | type<br>[big] trig pair<br>run5 run11 vent<br>vtac vfib emd<br>asys | -  | -   | -                   |

|  |  |                                |   |   |   |   |
|--|--|--------------------------------|---|---|---|---|
|  |  | category<br>[suprve]<br>condar | type<br>[1-bk] 2-b2<br>2-b2 3-bk<br>rbbb lbbb | - | - | - |
|--|--|--------------------------------|---|---|---|---|

**Note:**

The parameter in the lines that are adjacent to one another can be combined at will. For example, this allows for the function "sinus rhythm" alone more than 3,500 different options.

The abbreviations found under the main function "ECG" are listed and explained in alphabetical order in the following table:

| Display              | Meaning   |
|----------------------|---|
| arrhyth              | Arrhythmia  |
| condar               | Conduction Arrhythmia   |
| ECG form:<br>IEC 601 | Calibrating and analytical ECGs according to DIN EN 60601:2:51                            |
| ECG form:<br>others  | Other ECG forms (sinus rhythm ECGs and arrhythmia)  |
| HOR                  | horizontal heart position   |
| INT                  | intermediate heart position   |
| prem                 | premature arrhythmia  |
| QRS (mV)             | QRS Amplitude; for the norm ECGs calculated from zero, for all others as peak-peak values |
| QRS form             | QRS form, e.g., RS, R, QS   |
| sinusrh              | sinus rhythm  |
| STdev                | ST deviation  |
| suprve               | supraventricular arrhythmia   |
| type<br>analytical   | analytical ECGs, derivations are different  |
| type<br>calibration  | Calibrating ECGs; derivations are the same. Derivation III is a zero line.                |

|              |                         |
|--------------|-------------------------|
| <b>VERT</b>  | vertical heart position |
| <b>vntar</b> | ventricular arrhythmia  |

## List of Arrhythmia

| Display                 | Meaning                             |
|-------------------------|-------------------------------------|
| <b>category suprv</b>   |                                     |
| af1                     | atrial fibrillation, coarse         |
| af2                     | atrial fibrillation, fine           |
| aflt                    | atrial flutter                      |
| mb120                   | missed beat at 120 bpm              |
| mb80                    | missed beat at 80 bpm               |
| nodal                   | nodal rhythm                        |
| pat                     | paroxymal atrial tachycardia        |
| sa                      | sinus arrhythmia                    |
| svt                     | supraventricular tachycardia        |
| <b>category prem</b>    |                                     |
| mf1-3                   | multifocal PVCs                     |
| pac                     | premature atrial contraction        |
| pnc                     | premature nodal contraction         |
| pvc                     | premature ventricular contraction   |
| <b>category [vntar]</b> |                                     |
| asys                    | asystole                            |
| big                     | bigeminy                            |
| emd                     | electromotive disassociation        |
| pair                    | couplet (2 PVCs)                    |
| run11                   | 11 PVCs, 8 normal beats             |
| run5                    | 5 PVCs, 8 normal beats              |
| trig                    | trigeminy                           |
| vent                    | ventricular rhythm                  |
| vfib                    | ventricular fibrillation            |
| vtac                    | ventricular tachycardia             |
| <b>category suprv</b>   |                                     |
| 1-bk                    | First degree A-V block              |
| 2-b2                    | Second degree A-V block, Wenckebach |
| 2-b2                    | Second degree A-V block             |
| 3-bk                    | Third degree A-V block              |
| lbbb                    | left bundle branch block            |
| rbbb                    | right bundle branch block           |

## Main Function "RESP" - Respiration

In the basic setting after turning on the PSI 300, after pressing the main function key "RESP," the following settings in Lines 2 to 4 are shown in the display:

|        |        |     |      |      |           |
|--------|--------|-----|------|------|-----------|
| ECG    | [RESP] | BP  | TEMP | C.O. | NS        |
| Wave   | [Norm] |     |      |      |           |
| Rate   | 15     | 20  | [30] | 40   | 60 120    |
| R delt | 0.0    | 0.1 | 0.2  | 0.5  | [1.0] 3.0 |

Line 1 (Main Menu)

Line 2

Line 3

Line 4



### Note:

The settings and parameter in Line 5 to Line 7 can be shown and changed by pressing the "down" cursor key.

The basic settings of the main function "ECG" can be altered according to the following table. Each setting can be comined with all other settings.

| Line 2      | Line 3                          | Line 4                                 | Line 5                            | Line 6         | Line 7                           |
|-------------|---------------------------------|--|-----------------------------------|----------------|----------------------------------|
| Wave [Norm] | Rate<br>15 20 [30] 40<br>60 120 | R delt<br>0.0 0.1 0.2 0.5<br>[1.0] 3.0 | Ratio<br>1/1 2/1 [3/1]<br>4/1 5/1 | Lead<br>I [II] | Basel<br>500 [1000]<br>1500 2000 |

The abbreviations used in the main function "RESP" are shown and explained in alphabetical order in the following table:

| Display | Meaning  |
|---------|--|
| Wave    | Wave progression   |
| Norm    | Normal (spontaneous breathing)                               |
| Rate    | Breathing frequency [pro minute]                             |
| R delt  | Delta resistance (breath amplitude, resistance change [Ohm]) |
| Ratio   | Inspiration/Expiration Ratio                                 |
| Lead    | ECG Connection with which respiration can be measured        |
| Basel   | Baseline (simulated basic resistance of the thorax [Ohm])    |

## Main Function "BP" – Invasive Blood Pressure

In the basic setting after turning on the PSI 300, after pressing the main function key "BP," the following settings in Lines 2 to 4 are shown in the display:

|      |        |      |           |      |    |                    |        |
|------|--------|------|-----------|------|----|--------------------|--------|
| ECG  | RESP   | [BP] | TEMP      | C.O. | NS | Line 1 (Main Menu) |        |
| Chan | [IBP1] | IBP2 | IBP3      | IBP4 |    | Line 2             |        |
| Type | [Wave] | Stat | Swan-Ganz |      |    | Line 3             |        |
| Wave | ART    | LV   | [CVP]     | RV   | PA | PAW                | Line 4 |

### Note:

Further lines in this main function **cannot** be shown and changed by pressing the "down" cursor key.

The basic settings of the main function "BP" can be altered according to the following table.

| Line 2                                  | Line 3             | Line 4   | Note  |
|---|--------------------|--|---|
| Chan:<br>[IBP1]<br>IBP2<br>IBP3<br>IBP4 | Type:<br>[Wave]    | Wave:<br>ART LV [CVP] RV PA PAW                      |   |
|   | Type:<br>Stat      | Stat:<br>-10 -5 0 20 40 [80] 100                     |   |
|   | Type:<br>Swan-Ganz | Current Wave:<br>[ATM] insert remove                 | Selection with Insert:<br>⇒ CVP ⇒ RV ⇒ PA ⇒ PAW,<br>Selection with Remove:<br>Back to [ATM] |
|   |                    | Current Wave:<br>[PAW] deflate ⇔ [PA] inflate remove | Selection with Deflate:<br>Tp [PA] and selection: inflate /<br>remove                       |

The abbreviations used in the main function "BP" are shown and explained in alphabetical order in the following table:

| Display                   | Meaning   |
|---------------------------|---|
| <b>Chan:</b>              | Blood pressure channel for the current choice of setting          |
| Chan:<br><b>IBP1..4</b>   | Channels 1 to 4 (Jacks BP1 to BP4)                                |
| <b>CuWa:</b>              | Current Wave (current curve of the Swan-Ganz catheter function)   |
| CuWa:<br><b>ATM</b>       | Atmospheric Pressure (zero line, catheter outside of the patient) |
| CuWa:<br><b>CVP</b>       | Central Venous Pressure   |
| CuWa:<br><b>RV</b>        | Right Ventricular Pressure  |
| CuWa:<br><b>PA</b>        | Pulmonary Artery Pressure   |
| CuWa:<br><b>PAW</b>       | Pulmonary Artery Wedge Pressure                                   |
| <b>insert</b>             | Swan-Ganz Catheter, continue to insert                            |
| <b>remove</b>             | Swan-Ganz Catheter, remove  |
| <b>inflate</b>            | Swan-Ganz Balloon, inflate  |
| <b>deflate</b>            | Swan-Ganz Balloon, deflate  |
| <b>Type</b>               | Type of Blood Pressure  |
| Type:<br><b>Stat</b>      | Static (static pressure)  |
| Type:<br><b>Wave</b>      | Dynamic Pressure Courses  |
| Type:<br><b>Swan-Ganz</b> | Swan-Ganz Procedure   |
| Wave:<br><b>ART</b>       | Arterial Pressure   |
| Wave:<br><b>ATM</b>       | Atmospheric Pressure (zero line)                                  |
| Wave:<br><b>CVP</b>       | Central Venous Pressure   |
| Wave:<br><b>LV</b>        | Left Ventricular Pressure   |
| Wave:<br><b>PA</b>        | Pulmonary Artery Pressure   |
| Wave:<br><b>PAW</b>       | Pulmonary Artery Wedge Pressure                                   |
| Wave:<br><b>RV</b>        | Right Ventricular Presssure                                       |

## Main Function "TEMP" - Temperature

In the basic setting after turning on the PSI 300, after pressing the main function key "TEMP," the following settings in Lines 2 and 3 are shown in the display:

|     |      |      |        |      |      |       |    |    |     |
|-----|------|------|--------|------|------|-------|----|----|-----|
| ECG | RESP | BP   | [TEMP] | C.O. | NS   |       |    |    |     |
| FIX | 34   | 35   | 36     | [37] | 38   | 39    | 40 | 41 | 42> |
| VAR | 34   | [37] | 40     | Hypo | Hypr | Spike |    |    |     |
|     |      |      |        |      |      |       |    |    |     |

Line 1 (Main Menu)

Line 2

Line 3

Line 4

### Note:

Further lines in this main function **cannot** be shown and changed by pressing the "down" cursor key.

The basic settings of the main function "TEMP" can be altered according to the following table:

- **FIX (Channel Temp1)**

Channel Temp 1 should be used for more exact measurements because it is more precise (deviation of 0.1K).

- **VAR (Channel Temp2)**

Channel Temp2 should be used for sequences with lower precision (deviation of 0.4K):

- **Hypo** (Hypotrophy change of temperature)

Decreases temperature from 37°C to 34°C in 20 minutes, holds constant at 34°C for 60 minutes, then returns to 37°C in 20 minutes, holding at 37°C for 60 minutes

Continuous repetition of this interval.

- **Hypr** (Hypertrophy change of temperature)

Increases temperature from 37°C to 40°C in 20 minutes, holds at 40°C for 60 minutes, then decreases temperature from 40°C to 37°C in 20 minutes, holding at 37°C for 60 minutes.

Continuous repetition of this interval.

- **Spike** (Spike temperature changes)

Increases temperature from 37°C to 40°C in 20 minutes, then decreases back to 37°C in 20 minutes, holding at 37°C for 60 minutes.

Continuous repetition of this interval.



| Line 2                                    | Line 3                            | Line 4 | Line 5 |
|---|-----------------------------------|--------|--------|
| FIX<br>34 35 36 [37] 38 39 40 41<br>42 43 | VAR<br>34 [37] 40 Hypo Hypr Spike | -      | -      |

The abbreviations used in the main function "TEMP" are shown and explained in alphabetical order in the following table:

| Display      | Meaning                           |
|--------------|-----------------------------------|
| <b>FIX</b>   | Adjust Channel Temp1              |
| <b>VAR</b>   | Adjust Channel Temp2              |
| <b>Hypo</b>  | Hypotrophy Change of Temperature  |
| <b>Hypr</b>  | Hypertrophy Change of Temperature |
| <b>Spike</b> | Spike Temperature Changes         |

## Main Function "C.O." – Cardiac Output

In the basic setting after turning on the PSI 300, after pressing the main function key "CO," the following settings in Lines 2 and 3 are shown in the display:

|                 |      |     |      |        |           |
|-----------------|------|-----|------|--------|-----------|
| ECG             | RESP | BP  | TEMP | [C.O.] | NS        |
| Blood temp (°C) |      | 36  | [37] | 38     |           |
| min vol (l/min) |      | 3.0 | 4.0  | 4.5    | [5.0] 5.> |

Line 1 (Main Menu)  
 Line 2  
 Line 3  
 Line 4

### Note:

Further lines in this main function **cannot** be shown and changed by pressing the "down" cursor key.

The basic settings of the main function "CO" can be altered according to the following table:

| Line 2                       | Line 3  | Line 4 | Line 5 |
|------------------------------|---|--------|--------|
| Blood temp(°C)<br>36 [37] 38 | min vol (l/min)<br>3.0 4.0 4.5 [5.0] 5.5<br>6.0 7.0 intr lvs slow | -      | -      |

The abbreviations used in the main function "CO" are shown and explained in alphabetical order in the following table:

| Display                         | Meaning                           |
|---------------------------------|-----------------------------------|
| Blood temp (°C)                 | Blood Temperature                 |
| min vol (l/min):                | Cardiac output per minute [l/min] |
| min vol (l/min):<br><b>intr</b> | Interrupted Injectate Waveform    |
| min vol (l/min):<br><b>lvs</b>  | Left Ventricular Shunt            |
| min vol (l/min):<br><b>slow</b> | Slow Injectate Waveform           |

## Main Function "NS" - Noise

In the basic setting after turning on the PSI 300, after pressing the main function key "NS," the following settings in Line 2 are shown in the display:

|       |       |      |      |      |          |
|-------|-------|------|------|------|----------|
| ECG   | RESP  | BP   | TEMP | C.O. | [NS]     |
| Noise | [off] | 50Hz | 60Hz | HF   | Bas B+HF |
|       |       |      |      |      |          |
|       |       |      |      |      |          |

Line 1 (Main Menu)  
 Line 2  
 Line 3  
 Line 4

### Note:

Further lines of the main function "NS" **cannot** be shown and changed by pressing the "down" cursor key.

The basic settings of the main function "NS" can be altered according to the following table:

| Line 2                                | Line 3 | Line 4 |
|---------------------------------------|--------|--------|
| Noise:<br>[off] 50Hz 60Hz HF Bas B+HF | -      | -      |

The various disturbances are those that are inherent in DIN EN 60601-2-51. Thus, the filters of the ECG monitors can be verified.

The abbreviations used in the main function "NS" are shown and explained in alphabetical order in the following table:

| Display | Meaning   |
|---------|---|
| off     | No interferences                                    |
| 50Hz    | 50 Hz disturbance on the ECG                        |
| 60Hz    | 60 Hz disturbance on the ECG                        |
| HF      | High frequency disturbances                         |
| Bas     | Baseline digression (Frequency: 0.3 Hz)             |
| B+HF    | Baseline digression and high frequency disturbances |

### 13. Environmental Factors

Please always adhere to the operating and storage conditions of the PSI 300 in order to ensure that it functions without flaws:

- **Operating Conditions**

The PSI 300 should be operated only between the temperatures +5°C and +40°C and at a relative humidity between 50% and 80% (noncondensing).

- **Storage Conditions**

Store the PSI 300 only at temperatures between -20°C and +70°C at a relative humidity between 35% and 90% (noncondensing).

*When not in use over a longer period of time, remove the batteries.*

### 14. Cleaning and Disinfection

- **Cleaning**

The surfaces can be cleaned with any mild commercial plastic cleanser and a damp cloth.

**Caution!**

Do not allow any liquids to leak through the openings or jacks to the inside of the apparatus.

- **Desinfection**

Use any commercial disinfectants.

**Caution!**

Before using any cleansers and disinfectants, always test that they do not damage the plastic surfaces.

**Caution!**

***The PSI 300 is not suited to sterilization!***

## 15. Errors, Defects and Customer Service

### Error Codes

The PSI 300 is equipped with internal monitoring of the timing for the microcontroller and an internal voltage test. If the system detects an error, an error code appears on the display:

- Code 22 – "Timing Error"
- Code 25 – "Error ECG Voltage"

**If the error code is still displayed after a cold start, there is presumably another defect in the apparatus. In this case, please contact mtk biomed.**

### Trouble-Shooting

- **The PSI 300 cannot be switched on.**
  - Check whether the batteries are inserted correctly (position and polarity).
  - Check whether the batteries are sufficiently loaded.
- **The PSI 300 shuts itself off after a short time.**
  - Check whether the batteries are sufficiently loaded.
- **It's hard to read the display of the PSI 300.**
  - Check whether the batteries are sufficiently loaded.
- **The display shows a menu segment, but the keys don't react.**
  - An internal security shutdown has taken place. Switch off the PSI 300 for about 60 seconds and then switch it on again.
- **The display shows the menu segments, but no signals reach the jacks.**
  - An internal security shutdown has taken place. Switch off the PSI 300 for about 60 seconds and then switch it on again.

**If none of the above are able to remedy the problem, presumably there is another defect in the apparatus. In this case, please contact mtk biomed.**

## 16. Maintenance

The PSI 300 is maintenance-free. The output signals and parameter are

regulated by a microcontroller in the operating software (firmware). There is no adjustment hardware for alignment inside the PSI 300.

### Updating the Operating Software (Firmware Update)

A firmware update can be necessary if you want to expand the system of your PSI 300. In this case, please contact mtk biomed about a firmware update.

### Determining the Current Software

The version and the firmware revision (xx) of the PSI 300 are shown in the display for about 2 seconds during the initialization process, immediately after switching it on.

### 17. Intervals for Calibration

The manufacturer recommends to calibrate the **PSI 300** every 12 months. Keep your calibration protocols handy.

### 18. Packaging Materials and Disposal

Dispose of the packaging materials according to your local laws and regulations. When the PSI 300 has reached the end of its lifespan, do not dispose of it with the regular household waste.



Clean and disinfect the PSI 300 before disposing of it. For information on environmentally safe disposal, contact your customer service.

## 19. Pictographs

For the meaning of the pictographs on the identification plate of the PSI 300, see the following table:

- "Observe the Instruction Manual"



*Before starting the PSI 300 the first time, read this instruction manual. Comply with the safety regulations.*

- "Not Household Waste"



*Do not dispose in the regular household waste.*

- "CE Symbol"



*The CE symbol of the manufacturer assures that the PSI 300 conforms to EU Regulations 73/23/EWG.*

## 20. Technical Specifications

**Measurements and Weight** ca. 22.5 cm x 22 cm x 9 cm, 1.1 kg (without batteries)

### Power

- Rated Voltage: Batteries: 4 x 1.5V or Storage Batteries: 4 x 1.2V  
- Operating Voltage: 3.8 V= ... 6.8 V=  
- Power Input: 28 mA .... 46 mA, 38 mA typically

**Interface:** Serial Interface RS-232,  
for external steering via PC

### Environmental Conditions

#### - In Use

Temperature, relative humidity: +0°C ... +40 C, 35% ... 80% noncondensing

#### - Transport and Storage

Temperature, relative humidity: -20°C ... +70°C, 35% ... 90% noncondensing

### Performance Features

#### ECG:

12-Channel ECG with individual jacks for electrodes (10 electrodes R, L, F, N, C1 to C6)  
Bipolar Extremity Output according to Einthoven: I, II, III;  
Unipolar Extremity Output according to Goldberger: aVR, aVL, aVF

Generator lead according to Wilson: V1 - V6

#### Respiration:

Impedance Pneumograph, Output via ECG – Electrodes (choice of L or F)

#### Invasive Blood Pressure:

4 Blood pressure channels, able to be programmed individually  
Pneumograph Sensitivity 5µV/V/mmHg, utilizing the monitor reference voltage

#### Temperature

Simulation of the standard temperature sensor YSI400 via 2 channels  
Fixed adjustment (Temp1) and temperature sequences (Temp2)

#### Cardiac output

1 channel, measurement according to the thermo dilution principle



## Output Signals and Parameter

### ECG

- Calibrating and analytical ECG: DIN EN 60601-2-51:2004
  
- Sinus rhythm ECGs:
  - Choice of 7 frequencies:  
30, 60, 80, 120, 160, 200, 240
  - 3 different types of positions:  
intermediate, horizontal, vertical
  - 12 different amplitudes:  
0.05, 0.1, 0.2, 0.5, 0.7, 1.0, 1.5, 2.0, 2.5, 3.0, 4.0, 5.0 mV
  - ST Elevation and Sinking in 13 variations:  
+/- 800, +/- 600, +/- 400, +/- 200, +/- 100, +/- 50, 0  $\mu$ V
  - All parameter can be chosen independently and can be combined freely with one another.
  - Additional neonatal ECG, with frequencies, type of position, and amplitudes freely interchangeable .
  
- Arrhythmia:
  - auricular fibrillation, auricular flutter, extrasystole, bigeminy, ventricular tachycardia,
  
  - AV blocks of varying degrees, right bundle branch block (Rechtsschenkelblock)
  
- Supraventricular arrhythmia: 9 varieties
- Premature beats: 6 varieties
- Ventricular arrhythmia: 10 varieties
- Impulse lead error: 6 varieties

### Respiration

- 6 various frequencies,
- 6 various amplitudes,
- 5 various inhale/exhale ratios
- All parameter can be chosen independently and can be combined freely with one another,
- choice of Output I or Output II lying down, with 4 different baselines, also independent choice

## Invasive Blood Pressure

- static pressures (per channel): 7 various pressures  
(-10, 5, 0, 20, 40, 80, 100 mmHg)
- dynamic blood pressure sequences: artery pressure (ART)  
left ventricular pressure (LVP),  
central venous Pressure (CVP),  
right ventricular pressure (RVP),  
pulmonary pressure (PAP),  
pulmonary artery closing pressure (PCWP)
- Swan-Ganz Procedure: Swan-Ganz Procedure with different stations

## Temperature

- Starting temperatures: 34°C to 43 °C, adjustable in 1K steps
- Continual temperature sequences: hypertroph or hypotroph, on different channels

## Cardiac Output

- Flow sequence curve: after choosing the blood reference temperature

## Precision

- ECG: Frequency  $\pm 1\%$ , Amplitude  $\pm 5\%$
  - Respiration: Frequency  $\pm 1\%$ , Amplitude  $\pm 10\%$
  - Invasive blood pressure: Frequency  $\pm 1\%$ , Amplitude  $\pm 5\%$
  - Temperature: Temp1 (1K-Schritte):  $\pm 0.1K$ , Temp2:  $\pm 0.4K$ , Fr:  $\pm 1\%$
  - Cardiac output: Frequency  $\pm 1\%$ , Amplitude  $\pm 5\%$
-